





Metallic Magnetic Calorimeters: Detectors for High Resolution X-Ray Spectroscopy




A. Fleischmann,
T. Daniyarov, M. Linck
H. Rotzinger, A. Burck, C. Enss
Kirchhoff-Institut für Physik
Universität Heidelberg

G.M. Seidel
Department of Physics
Brown University

Physics of MMCs
Recent Results
Applications

Metallic Magnetic Calorimeter



$$\delta M = \frac{\partial M}{\partial T} \delta T = \frac{\partial M}{\partial T} \frac{E_\gamma}{C_{\text{tot}}}$$

Sensor: Au:Er ←


Ag:Er

Bi₂Te₃:Er

PbTe:Er

Magnetism of Sensor Material: Au:Er

gold (fcc)



4f - electrons (Er)

Rare earth ions in metals


- 4f - shell partially filled




angular momentum & magnetic moment

Er^{3+} : $J = 15/2 \rightarrow 16$ levels

- crystal field splitting




Magnetic Susceptibility of Au:Er



Calculation of thermodynamic properties


Interactions:

- magnetic dipole-dipole
- RKKY
(indirect exchange interaction)





Numerical simulation:


- 12 spins, randomly distributed
- Hamiltonian matrix
 - ↓ numeric diagonalization
- energy eigenvalues
 - ↓ $M(B,T)$, $C(B,T)$, 10.000 configs
- $M(B,T)$, $C(B,T)$



Magnetization and Heat Capacity





$$\delta \Phi_S \propto \frac{\partial M}{\partial T} \frac{1}{C_{\text{tot}}} \delta E$$



Calorimeter Signal

$$\delta\Phi_S = f(r, h) \frac{\partial M}{\partial T} \frac{1}{C_{\text{tot}}} \delta E$$



- satisfying agreement of theory and experiment
- signal size can be predicted!
- optimal detector design?

$$S_{\text{max}}(C_a, g, \alpha, T, x, r, h, B)$$

Optimized Detector Design


- Application requires:
 - Absorber size \rightarrow heat capacity
 - Operating temperature
- Sensor material (e.g. Au:Er):
 - Landé-factor: g
 - RKKY-interaction: a

• Optimization of sensitivity

\rightarrow Optimal sensor parameters:

- diameter, thickness
- spin concentration
- magnetic field

Optimized Sensor:



- Sensitivity ($S = \Delta\Phi/\Delta E$):


$$S_{\max} = 1,2 \times 10^{-9} \left(\frac{g}{C_{\text{Abs}} \propto T^2} \right)^{1/3} \frac{\Phi_0}{\text{eV}}$$

$$S_{\max} \propto \frac{1}{T^{2/3}} \quad \text{increases with decreasing } T$$

$$S_{\max} \propto \frac{1}{C_{\text{Abs}}^{1/3}} \quad \text{weakly dependent on absorber heatcapacity}$$


Alternative Geometries

coupling depends on spin position



$$S_{\max} \propto \frac{1}{C_{\text{Abs}}^{1/3}}$$

Meanders, multiple loops, miniature coils, ...



visit Poster A 08

Sensitivity dF/dE :


- almost independent of absorber heat capacity
- roughly proportional to the inverse of structure width

2eV X-Ray Detector Pixel



- Temperature: 50mK
- Absorber: $250 \times 250 \times 5 \mu\text{m}^3$, gold
 - heat capacity: 1 pJ/K
 - quantum efficiency: $> 98\%$ (6keV)


- ➔ • Sensor:
 - Au:Er 900ppm, $\varnothing 52\mu\text{m}$, $h=6.5\mu\text{m}$
- Sensitivität:
 - $0,12 \Phi_0 / \text{keV}$
- Relaxation times:
 - rise: $< 10^{-6}\text{s}$
 - decay: 10^{-3}s



- ➔ • Energy resolution: $\Delta E_{\text{FWHM}} \cong 1.4 \text{ eV}$

SQUID noise
and Johnson noise:
20% increase

Two-Stage SQUID Magnetometer



- Detector SQUID read out by a **current-sensor** SQUID
- **Two-stage setup** allows for:
 - **Low noise**
 - **Large bandwidth** / Slewrate
 - **Small power dissipation** on detector SQUID chip (Voltage bias)



Two-Pixel Detector

SQUID Gradiometer With Two Sensors

two Au:Er 300 ppm sensors

\varnothing 50 μm , $h = 25 \mu\text{m}$

magnetic field 3 mT

two gold absorbers: $160 \times 160 \times 5 \mu\text{m}^3$

Detector SQUID:

commercial SQUID chip
M.B. Ketchen, IBM 1992


Amplifier SQUID: IPHT-Jena

SQUID-Electronics (IPHT-Jena):

voltage noise: $0.3 \text{ nV}/\sqrt{\text{Hz}}$
very low temperature dependence


ADR: VeriCold Technologies, Munich

21 mK base temperature, holdtime below 30 mK: 2 days




Two-Pixel Detector

X-rays from ^{55}Fe -Source:



➔ clean spectrum

Two-Pixel Detector





performance of pixels
almost identical

Resolution: K_{α} -Line ^{55}Mn

Two Au:Er 300 ppm sensors

Two Gold absorbers: $160 \times 160 \times 5 \mu\text{m}^3$

Total heatcapacity: 2.5 pJ/K



energy resolution 3.4 eV



resolving power > 1700

Resolution: K_{β} -Line ^{55}Mn

Two Au:Er 300 ppm sensors


Two Gold absorbers: $160 \times 160 \times 5 \mu\text{m}^3$

Total heatcapacity: 2.5 pJ/K




energy resolution 3.4 eV

Baseline Noise




Resolution is ,constant'
over the whole energy range
Small temperature drifts
Small position dependence

Noise Spectrum at 33 mK



Noise at low frequencies:
about 2 times larger
than expected from
thermodynamic energy fluctuations

1/f - Noise



Observations:

temperature **independent**
(30 mK ... 4 K)

magnetic field **independent**
(0 ... 35 G)

concentration dependence ?

volume dependence ?


Contribution to the energy resolution:
2 eV



Thermalization


Sensor bonded with vacuum grease to Si

Heat capacity 1.2×10^{-12} J/K



Sensor: 1000x larger,
spot-welded to copper block

Heat capacity 10^{-9} J/K



- magnetic calorimeters can be made very fast
- Plan: microstructured thin film on the SQUID-chip for thermalization

Applications

Industrial :

material analysis
absolute activity measurement

Astrophysics:

solar neutrino detection
x-ray astronomy
dark matter detection


Nuclear physics:

gamma-ray detection
beta decay
double-beta decay



Material Analysis

X-ray fluorescence analysis



long term goal is: resolution $2 \text{ eV at } 6 \text{ keV}$,
count rate $100/\text{s}$,
quantum efficiency $99\% \text{ up to } 6 \text{ keV}$
operated in a **cryo-liquid free cryostat**



MMCs for Metrology

Goal:

Measurement of absolute source activity

Challenge:

99.99% quantum efficiency

4 π -geometry

large dynamic range



details on Poster N
04, (M. Loidl, et al.)




COMMISSARIAT À L'ÉNERGIE ATOMIQUE

Example ^{55}Fe

4 π -Sandwich:



Recent Result




Low energy background


→ shielding & higher activity



Summary



- MMCs can be used for
high resolution x-ray spectroscopy
- MMCs are a promising detector concept
especially for applications which
require high heat capacities



- There are lots of things to be looked at
- Hopefully more groups will join this effort ... ?

